

Rinofast *aerosol*



- Nasal congestion
- Allergic and non-allergic rhinitis
- Rhinosinusitis
- Mucociliary clearance
- Nasal hygiene

Medical Device CE



MADE IN ITALY

NASAL CONGESTION & COLD

PLUS

DESCRIPTION

Rinofast aerosol is a hypertonic saline solution (3% sodium chloride) at pH controlled, ready to use. Thanks to its osmotic action, Rinofast Aerosol promotes the restoration of normal mucociliary clearance, allows fluidization, transport, and removal of mucus, favors nasal decongestion improving respiratory function. It helps in reduction of dryness gently hydrating nose and mouth. Suitable for children, also in case of colds.



NASAL IRRIGATION: half vial in each nostril 2-3 times a day
INHALATION THERAPY: 1 vial, 1-2 times a day

THERAPEUTIC INDICATIONS

- Cold
- Nasal congestion
- Nasal hygiene

PRESENTATION

Each box of Rinofast Aerosol contains 15 resealable vials of 5 ml.

INGREDIENTS

Water, Sodium Chloride 3%, Potassium Citrate, Citric Acid.

REGULATORY STATUS

- Class II
- Free sale certificate
- CE-Conformity certificate
- Certificate of origin (at the time of invoicing)
- ISO 13485 certificate



DOUBLE ADMINISTRATION ROUTE





NEBULIZER



DIRECT IRRIGATION

NaCl
Immediate action

- Osmotic action
- Improves mucociliary clearance
- Removes the mucus



reusable within
12h
reusable within

RESEALABLE VIALS



CLINICAL TEST

- **RISK ANALYSIS** - according to UNI CEI EN ISO 14971
- **BIOCOMPATIBILITY** - according to UNI EN ISO 10993: CYTOTOXICITY for DIRECT CONTACT - according to UNI EN ISO 10993-5:2009
- **NASAL MUCOSA IRRITATION TEST** - according to UNI EN ISO 10993-10:2010
- **DELAYED HYPERSENSIVITY TEST (GPMT)** - according to UNI EN ISO 10993-10:2010
- **CLINICAL EVALUATION**
- **STABILITY** - according to ICH Guidelines: ACCELERATED STABILITY STUDY
AFTER OPENING ANALYSIS-IN USE TEST

Medical Device CE

